

EVH System Checklist

Equipment	
	Video Tower (monitor, camera, light source)
	CO ₂ Source – ie. Tank – with filter and regulator*
	Sterile CO ₂ tubing
	Sterile camera
	Sterile light cord
	Sterile 5mm, 30° Scope (29cm in length)
	Bipolar Generator with Foot Pedal
	Sterile Bipolar Cord
	ClearGlide EVH Kit
	Endoloop/TTS Loop
Optional	
	Marking Pen
	Sterile Suction Tip/Tubing
	Anti-Fog

* CO₂ Settings – approx. 6L/min
 Harmonic Scalpel Settings – Min = 3, Max = 5 – Refer to Harmonic Scalpel Handbook

ERA System Checklist

Equipment	
	Video Tower (monitor, camera, light source)
	CO ₂ Source – ie. Tank – with filter and regulator*
	Sterile CO ₂ tubing
	Sterile camera
	Sterile light cord
	Sterile 5mm, 30° Scope (29cm in length)
	Harmonic Scalpel with Foot Pedal
	Sterile Wrench
	ClearGlide ERA Kit
	Endoloop/TTS Loop
Optional	
	Marking Pen
	Sterile Suction Tip/Tubing
	Anti-Fog

* CO₂ Settings – approx. 6L/min

Harmonic Scalpel Settings – Min = 3, Max = 5 – Refer to Harmonic Scalpel Handbook